



Instrumentation Laboratory ILab 300 instructions

rev 426.0.2 - 2010-08-01

Chema
D I A G N O S T I C A

Chema Diagnostica
Via Campania 2/4
I-60030 Monsano | An Italy
T +39 0731 605064 | F +39 0731 605672
www.chema.com | mail@chema.com

Instrumentation Laboratory ILab 300 instructions

rev 426.0.2 - 2010-08-01

Warning!

These programs should be used only as guideline. To check the reagent correct setting on the instrument, use only good quality control sera and work accordingly to with good laboratory practice. Instrumental factor must be checked. Chema Diagnostica will not be responsible for instrumental bad programming.

Analyzer: ILab 300

Application: **ACID PHOSPHATASE (TOTAL)** - Codes AC 0120 TC
 Preparation: POWDER SINGLE REAGENT.
 FOLLOW PREPARATION PROGRAM INDICATED IN INSERT SHEET.
 Storage: REFRIGERATE AT 2-8°C
 Stability: 10 DAYS

PROGRAM

Description: ACP
 Unit: UI
 Decimals: 1
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	0.0	0.0	0.0	0.0	4.7	75.0	75.0	75.0
Female:								
Children:								
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No				No	No	

Reaction type: Kinetic
 Direction: Up
 E.P. limit (abs): N/A
 Depl.limit (abs): 2.0000
 First limit (abs): N/A
 Linear factor: N/A
 Fit: 0.950

		Parameter					
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):			0	0	0	242	252
Dil./Rgt. Code:			ACP	0			
Lot Number:							
Ratio/Vol (µl):	1/1		300	0	0		
Rinse (µl):			0	0	0		
Sample (µl):			30				

Filter 1 (nm): 405
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 0.0 RBL min (abs): -0.2000
 High: 75.0 Max (abs): 2.0000

Calculation. model: Factor RBL stability (days): 1 ** definito dall'operatore
 Factor: 743.00 Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: NESSUNO
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **ALBUMIN** - Codes BC 0100/ 0500 / 1000 / 1500 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: ROOM TEMPERATURE (2-30°C)
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

Description: ALBUMIN
 Unit: g/dl
 Decimals: 1
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

Reference Range

	LOW VALUES				HIGH VALUES			
Male:	1.0	1.0	1.0	3.8	5.1	7.0	7.0	7.0
Female:	1.0	1.0	1.0	3.8	5.1	7.0	7.0	7.0
Children:	1.0	1.0	1.0	3.8	5.1	7.0	7.0	7.0
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Parameter

Reaction type: End Point	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read	
Direction: Up		0	0	0	98	0	
E.P. limit (abs): 0.1000		ALB	0				
Depl.limit (abs): N/A		Lot Number:			Filter 1 (nm): 620 Filter 2 (nm): (none) Bicr. factor: 1.00		
First limit (abs): N/A	1/1	Ratio/Vol (µl):	400	0			0
Linear factor: N/A		Rinse (µl):	0	0			0
Fit: N/A		Sample (µl):	3				

Lin limit. low: 0.2 RBL min (abs): -0.2000
 High: 7.0 Max (abs): 0.8000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 7 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: AUTOCAL H
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **ALKALINE PHOSPHATASE FL (DGKC)** - Codes AL F080 / F245 / F400 / F600 CH
 Preparation: MIX REAGENTS ACCORDING TO INSERT SHEET

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: ALK.PHOS.
 Unit: UI
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

Reference Range

	LOW VALUES				HIGH VALUES			
Male:	0	0	0	0	270	2500	2500	2500
Female:	0	0	0	0	240	2500	2500	2500
Children:	0	0	0	0	270	2500	2500	2500
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No				No	No	No

Parameter

	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Reaction type: Kinetic						
Direction: Up		0	0	0	44	144
E.P. limit (abs): N/A		ALP	0			
DepL.limit (abs): 2.0000		Lot Number:				
First limit (abs): N/A	1/1	Ratio/Vol (µl):	300	0		
Linear factor: N/A		Rinse (µl):	0	0		
Fit: 0.970		Sample (µl):	6			
					Filter 1 (nm): 405 Filter 2 (nm): (none) Bicr. factor: 1.00	

Lin limit. low: 0 RBL min (abs): -0.2000
 High: 2800 Max (abs): 1.5000

Calculation. model: Factor RBL stability (days): 1 ** definito dall'operatore
 Factor: 2757.00 Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: NESSUNO
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **ALKALINE PHOSPHATASE FL (IFCC)** - Codes AF F080 / F245 / F400 / F600 CH
 Preparation: MIX REAGENTS ACCORDING TO INSERT SHEET

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: ALK.PHOS.
 Unit: UI
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

Reference Range

	LOW VALUES			HIGH VALUES				
Male:	0	0	0	35	104	2500	2500	2500
Female:	0	0	0	40	129	2500	2500	2500
Children:	0	0	0	0	0	2500	2500	2500
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Parameter

	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Reaction type: Kinetic						
Direction: Up						
E.P. limit (abs): N/A						
DepL.limit (abs): 2.0000						
First limit (abs): N/A						
Linear factor: N/A						
Fit: 0.970						
Times (sec):		0	0	0	44	144
Dil./Rgt. Code:		ALP	0			
Lot Number:						
Ratio/Vol (µl): 1/1		300	0	0		
Rinse (µl):		0	0	0		
Sample (µl):		6				

Filter 1 (nm): 405
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 0 RBL min (abs): -0.2000
 High: 2800 Max (abs): 1.5000

Calculation. model: Factor RBL stability (days): 1 ** definito dall'operatore
 Factor: 2757.00 Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: NESSUNO
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **AMYLASE FL** - Codes AM F060 / F120 / F245 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: AMYLASE

Unit: UI
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

Reference Range

	LOW VALUES				HIGH VALUES			
Male:	0	0	0	0	96	2000	2000	2000
Female:	0	0	0	0	96	2000	2000	2000
Children:	0	0	0	0	96	2000	2000	2000
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Parameter

	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Reaction type: Kinetic		0	0	0	98	0
Direction: Up		AMYLASE	0			
E.P. limit (abs): N/A						
DepL.limit (abs): 2.0000						
First limit (abs): N/A	1/1	300	0	0		
Linear factor: N/A		0	0	0		
Fit: 0.970		6				
					Filter 1 (nm): 405 Filter 2 (nm): (none) Bicr. factor: 1.00	

Lin limit. low: 0 RBL min (abs): -0.2000
 High: 2000 Max (abs): 0.8000

Calculation. model: Factor RBL stability (days): 1 ** definito dall'operatore
 Factor: 3953.00 Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: NESSUNO
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **AMYLASE EPS FL** - Codes EA F080 / F245 CH
 Preparation: MIX REAGENTS ACCORDING TO INSERT SHEET

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: AMYLASE
 Unit: UI
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	0	0	0	28	100	1500	1500	1500
Female:	0	0	0	28	100	1500	1500	1500
Children:	0	0	0	28	100	1500	1500	1500
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Reaction type: Kinetic
 Direction: Up
 E.P. limit (abs): N/A
 Depl.limit (abs): 2.0000
 First limit (abs): N/A
 Linear factor: N/A
 Fit: 0.970

		Parameter					
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):			0	0	0	98	0
Dil./Rgt. Code:			AMYLASE	0			
Lot Number:							
Ratio/Vol (µl):	1/1		300	0	0		
Rinse (µl):			0	0	0		
Sample (µl):			10				

Filter 1 (nm): 405
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 0 RBL min (abs): -0.2000
 High: 1500 Max (abs): 0.8000

Calculation. model: Factor RBL stability (days): 1 ** definito dall'operatore
 Factor: 3480.00 Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: NESSUNO
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **PANCREATIC ISOAMYLASE FL** - Code PA F080 / F245 CH

Preparation: REAGENT # 1 - USE REAGENT A

REAGENT # 2 - USE REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: UNTIL EXPIRATION DATE ON LABEL

Description: ISOAMYLASE

Unit: UI

Decimals: 0

LIS Code: **

Unit Factor: 1.0

Slope: 1.00

Intercept: 0.00

Reference Range

	LOW VALUES			HIGH VALUES				
Male:	0	0	0	13	53	2500	2500	2500
Female:	0	0	0	13	53	2500	2500	2500
Children:	0	0	0	13	53	2500	2500	2500
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Parameter

	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Reaction type: Kinetic		0	206	0	62	170
Direction: Up		ISOAMY1	ISOAMY2			
E.P. limit (abs): N/A						
DepL.limit (abs): 2.0000						
First limit (abs): N/A	1/1	200	50	0		
Linear factor: N/A		0	0	0		
Fit: 0.970		5				
					Filter 1 (nm): 405 Filter 2 (nm): (none) Bicr. factor: 1.00	

Lin limit. low: 0

High: 2500

RBL min (abs): -0.2000

Max (abs): 0.8000

Calculation. model: Factor

Factor: 6280.00

Sample blank: No

RBL stability (days): 1

Calibration stab. (days): 99

Dynamic controls (min): None

** definito dall'operatore

N/A non applicabile

Calibratore suggerito: NESSUNO

Controlli suggeriti: QUANTINORM CHEMA

QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **BICARBONATE FL** - Code BR F060 / F245 / F400 CH
 Preparation: LIQUID REAGENT READY TO USE
 Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: BICARBONATE

Unit: mmol/L
 Decimals: 1
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

Reference Range

	LOW VALUES				HIGH VALUES			
Male:	2	2	2	22	29	50	50	50
Female:	2	2	2	22	29	50	50	50
Children:	2	2	2	20	28	50	50	50
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No				No	No	

Parameter

Reaction type: End Point
 Direction: Up
 E.P. limit (abs): 0.1000
 Depl.limit (abs): N/A
 First limit (abs): N/A
 Linear factor: N/A
 Fit: N/A

	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):		0	0	0	150	0
Dil./Rgt. Code:		CO2	0			
Lot Number:						
Ratio/Vol (µl):	1/1	360	0	0		
Rinse (µl):		0	0	0		
Sample (µl):		4				

Filter 1 (nm): 405
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 5 RBL min (abs): 1.0000
 High: 50 Max (abs): 2.5000

Calculation. model: Standard
 Factor: N/A
 Sample blank: No

RBL stability (days): 1
 Calibration stab. (days): 10
 Dynamic controls (min): None

** definito dall'operatore
 N/A non applicabile

Analyzer: ILab 300

Application: **TOTAL BILIRUBIN FL** - Code DT F125 / F500 CH

Preparation: REAGENT # 1 - USE REAGENT A
REAGENT # 2 - USE REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: UNTIL EXPIRATION DATE ON LABEL

Description: BILI-T

Unit: mg/dl

Decimals: 2

LIS Code: **

Unit Factor: 1.0

Slope: 1.00

Intercept: 0.00

Reference Range

	LOW VALUES			HIGH VALUES				
Male:	0.00	0.00	0.00	0.10	1.20	20.00	20.00	20.00
Female:	0.00	0.00	0.00	0.10	1.20	20.00	20.00	20.00
Children:	0.00	0.00	0.00	0.10	1.20	20.00	20.00	20.00
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Parameter

Reaction type: End Point

Direction: Up

E.P. limit (abs): 0.1000

Depl.limit (abs): N/A

First limit (abs): N/A

Linear factor: N/A

Fit: N/A

Times (sec):

Dil./Rgt. Code:

Lot Number:

Ratio/Vol (µl):

Rinse (µl):

Sample (µl):

Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
	0	180	0	440	0
	BT1	BT/D2			
1/1	300	75	0	<div style="border: 1px solid black; padding: 5px;"> Filter 1 (nm): 546 Filter 2 (nm): (none) Bicr. factor: 1.00 </div>	
	0	0	0		
	25				

Lin limit. low: 0.00

High: 20.00

RBL min (abs): -0.2000

Max (abs): 0.8000

Calculation. model: Standard

Factor: N/A

Sample blank: Yes

RBL stability (days): 1

Calibration stab. (days): 1

Dynamic controls (min): None

** definito dall'operatore

N/A non applicabile

Calibratore suggerito: AUTOCAL H

Controlli suggeriti: QUANTINORM CHEMA

QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **DIRECT BILIRUBIN FL** - Code DD F125 / F500 CH

Preparation: REAGENT # 1 - USE REAGENT A

REAGENT # 2 - USE REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: UNTIL EXPIRATION DATE ON LABEL

Description: BILI-D

Unit: mg/dl

Decimals: 2

LIS Code: **

Unit Factor: 1.0

Slope: 1.00

Intercept: 0.00

Reference Range

	LOW VALUES				HIGH VALUES			
Male:	0.00	0.00	0.00	0.00 0.20	20.00	20.00	20.00	20.00
Female:	0.00	0.00	0.00	0.00 0.20	20.00	20.00	20.00	20.00
Children:	0.00	0.00	0.00	0.00 0.20	20.00	20.00	20.00	20.00
Low Alert:	Low alert	Very low	Low	Normal values	High	Very high	High alert	
Rerun:	No	No				No	No	

Parameter

Reaction type: End Point

Direction: Up

E.P. limit (abs): 0.1000

Depl.limit (abs): N/A

First limit (abs): N/A

Linear factor: N/A

Fit: N/A

Times (sec):

Dil./Rgt. Code:

Lot Number:

Ratio/Vol (µl):

Rinse (µl):

Sample (µl):

Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
	0	180	0	62	0
	BD1	BT/D2			
1/1	300	75	0		
	0	0	0		
	25				

Filter 1 (nm): 546
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 0.00

High: 20.00

RBL min (abs): -0.2000

Max (abs): 0.8000

Calculation. model: Standard

Factor: N/A

Sample blank: Yes

RBL stability (days): 1

Calibration stab. (days): 1

Dynamic controls (min): None

** definito dall'operatore

N/A non applicabile

Calibratore suggerito: AUTOCAL H

Controlli suggeriti: QUANTINORM CHEMA

QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **CALCIUM** - Code CA 0305 / 0505 CH

Preparation: REAGENT # 1 - USE REAGENT B

REAGENT # 2 - USE REAGENT A

Storage: REFRIGERATE AT 2-8°C

Stability: UNTIL EXPIRATION DATE ON LABEL

Description: CALCIUM

Unit: mg/dl

Decimals: 2

LIS Code: **

Unit Factor: 1.0

Slope: 1.00

Intercept: 0.00

Reference Range

	LOW VALUES			HIGH VALUES				
Male:	1.00	1.00	1.00	8.80	10.80	20.00	20.00	20.00
Female:	1.00	1.00	1.00	8.80	10.80	20.00	20.00	20.00
Children:	1.00	1.00	1.00	10.00	11.60	20.00	20.00	20.00
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Parameter

Reaction type: End Point

Direction: Up

E.P. limit (abs): 0.1000

Depl.limit (abs): N/A

First limit (abs): N/A

Linear factor: N/A

Fit: N/A

Times (sec):

Dil./Rgt. Code:

Lot Number:

Ratio/Vol (µl):

Rinse (µl):

Sample (µl):

Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
	0	90	0	170	0
	CA1	CA2			
1/1	300	300	0	<div style="border: 1px solid black; padding: 5px;"> Filter 1 (nm): 578 Filter 2 (nm): (none) Bicr. factor: 1.00 </div>	
	0	0	0		
	15				

Lin limit. low: 0.80

High: 20.00

RBL min (abs): -0.2000

Max (abs): 2.0000

Calculation. model: Standard

Factor: N/A

Sample blank: Yes

RBL stability (days): 1

Calibration stab. (days): 1

Dynamic controls (min): None

** definito dall'operatore

N/A non applicabile

Calibratore suggerito: AUTOCAL H

Controlli suggeriti: QUANTINORM CHEMA

QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **CALCIUM ASX** - Code CA 0100 / 0500 CH
 Preparation: LIQUID REAGENT READY TO USE
 Storage: ROOM TEMPERATURE (2-30°C)
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: CALCIUM
 Unit: mg/dl
 Decimals: 2
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

Reference Range

	LOW VALUES				HIGH VALUES			
Male:	1.00	1.00	1.00	8.80	10.80	20.00	20.00	20.00
Female:	1.00	1.00	1.00	8.80	10.80	20.00	20.00	20.00
Children:	1.00	1.00	1.00	10.00	11.60	20.00	20.00	20.00
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No				No	No	

Parameter

	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Reaction type: End Point		0	0	0	170	0
Direction: Up		CALCIUM	0			
E.P. limit (abs): 0.1000						
Depl.limit (abs): N/A						
First limit (abs): N/A	1/1	400	0	0		
Linear factor: N/A		0	0	0		
Fit: N/A		4				
					Filter 1 (nm): 620 Filter 2 (nm): (none) Bicr. factor: 1.00	

Lin limit. low: 0.80 RBL min (abs): -0.2000
 High: 20.00 Max (abs): 2.0000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 1 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: AUTOCAL H
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **CHLORIDE** - Code CL 0100 / 0500 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: ROOM TEMPERATURE (2-30°C)
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: CHLORIDE
 Unit: meq/l
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	50	50	50	98	110	200	200	200
Female:	50	50	50	98	110	200	200	200
Children:	50	50	50	98	110	200	200	200
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Reaction type: End Point
 Direction: Up
 E.P. limit (abs): 0.1000
 Depl.limit (abs): N/A
 First limit (abs): N/A
 Linear factor: N/A
 Fit: N/A

		Parameter					
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):			0	0	0	368	0
Dil./Rgt. Code:			CL	0			
Lot Number:							
Ratio/Vol (µl):	1/5		450	0	0		
Rinse (µl):			0	0	0		
Sample (µl):			3				

Filter 1 (nm): 492
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 20 RBL min (abs): -0.2000
 High: 200 Max (abs): 0.8000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 7 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: AUTOCAL H
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **CHOLESTEROL FL** - Codes CT F100 / F400 / 100F / 150F CH

Preparation: LIQUID READY TO USE SINGLE REAGENT

Storage: REFRIGERATE AT 2-8°C

Stability: UNTIL EXPIRATION DATE ON LABEL

Description: CHOLESTEROL

Unit: mg/dl

Decimals: 0

LIS Code: **

Unit Factor: 1.0

Slope: 1.00

Intercept: 0.00

Reference Range

		LOW VALUES			HIGH VALUES			
Male:	5	5	5	140	200	700	700	700
Female:	5	5	5	140	200	700	700	700
Children:	5	5	5	140	200	700	700	700
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No				No	No	

Parameter

Reaction type: End Point

Direction: Up

E.P. limit (abs): 0.1000

Depl.limit (abs): N/A

First limit (abs): N/A

Linear factor: N/A

Fit: N/A

Times (sec):

Dil./Rgt. Code:

Lot Number:

Ratio/Vol (µl): 1/1

Rinse (µl):

Sample (µl):

Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
	0	0	0	350	0
	CHOL	0			
	360	0	0		
	0	0	0		
	4				

Filter 1 (nm): 510
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 5

High: 700

RBL min (abs): -0.2000

Max (abs): 0.8000

Calculation. model: Standard

Factor: N/A

Sample blank: No

RBL stability (days): 1

Calibration stab. (days): 10

Dynamic controls (min): None

** definito dall'operatore

N/A non applicabile

Calibratore suggerito: AUTOCAL H

Controlli suggeriti: QUANTINORM CHEMA

QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **HDL-DIRECT FL** - Code HD F080 / F245 / F400 CH

Preparation: REAGENT # 1 - USE REAGENT A

REAGENT # 2 - USE REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: UNTIL EXPIRATION DATE ON LABEL

Description: HDL-C

Unit: mg/dl

Decimals: 0

LIS Code: **

Unit Factor: 1.0

Slope: 1.00

Intercept: 0.00

Reference Range

	LOW VALUES			HIGH VALUES				
Male:	5	5	5	35	79	200	200	200
Female:	5	5	5	42	88	200	200	200
Children:	5	5	5	140	200	200	200	200
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Parameter

Reaction type: End Point

Direction: Up

E.P. limit (abs): 0.1000

Depl.limit (abs): N/A

First limit (abs): N/A

Linear factor: N/A

Fit: N/A

Times (sec):

Dil./Rgt. Code:

Lot Number:

Ratio/Vol (µl):

Rinse (µl):

Sample (µl):

Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
	0	252	0	242	0
	HDL1	HDL2			
1/1	270	90	0		
	0	0	0		
	3				

Filter 1 (nm): 620
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 2

High: 220

RBL min (abs): -0.2000

Max (abs): 2.0000

Calculation. model: Standard

Factor: N/A

Sample blank: Yes

RBL stability (days): 1

Calibration stab. (days): 7

Dynamic controls (min): None

** definito dall'operatore

N/A non applicabile

Calibratore suggerito: AUTOCAL H

Controlli suggeriti: QUANTINORM CHEMA

QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **LDL-DIRECT FL - Code DL F080 CH**

Preparation: REAGENT # 1 - USE REAGENT A

REAGENT # 2 - USE REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: UNTIL EXPIRATION DATE ON LABEL

Description: LDL-C

Unit: mg/dl

Decimals: 0

LIS Code: **

Unit Factor: 1.0

Slope: 1.00

Intercept: 0.00

Reference Range

	LOW VALUES			HIGH VALUES				
Male:	5	5	5	130	160	400	400	400
Female:	5	5	5	130	160	400	400	400
Children:	5	5	5	130	160	400	400	400
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Parameter

Reaction type: End Point

Direction: Up

E.P. limit (abs): 0.1000

Depl.limit (abs): N/A

First limit (abs): N/A

Linear factor: N/A

Fit: N/A

Times (sec):

Dil./Rgt. Code:

Lot Number:

Ratio/Vol (µl):

Rinse (µl):

Sample (µl):

Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
	0	252	0	242	0
	LDL1	LDL2			
1/1	270	90	0		
	0	0	0		
	3				

Filter 1 (nm): 620
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 2

High: 400

RBL min (abs): -0.2000

Max (abs): 2.0000

Calculation. model: Standard

Factor: N/A

Sample blank: Yes

RBL stability (days): 1

Calibration stab. (days): 7

Dynamic controls (min): None

** definito dall'operatore

N/A non applicabile

Calibratore suggerito: AUTOCAL H

Controlli suggeriti: QUANTINORM CHEMA

QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **CHOLINESTERASE FL (DGKC)** - Codes CH F096 / F245 CH

Preparation: REAGENT # 1 - USE REAGENT A

REAGENT # 2 - USE REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: UNTIL EXPIRATION DATE ON LABEL

Description: CHOLINESTERASE	Reference Range								
	LOW VALUES				HIGH VALUES				
Unit: UI	Male:	500	500	500	5600	11200	20000	20000	20000
Decimals: 0	Female:	500	500	500	4200	10800	20000	20000	20000
LIS Code: **	Children:	500	500	500	5600	11200	20000	20000	20000
Unit Factor: 1.0	Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Slope: 1.00	Rerun:	No	No				No	No	No
Intercept: 0.00									

Reaction type: Kinetic	Parameter					
	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Direction: Down	Times (sec):	0	180	0	126	98
E.P. limit (abs): N/A	Dil./Rgt. Code:	CHE1	CHE2			
DepL.limit (abs): 2.5000	Lot Number:					
First limit (abs): N/A	Ratio/Vol (µl):	1/1	300	60	0	
Linear factor: N/A	Rinse (µl):	0	0	0		
Fit: 0.950	Sample (µl):	6				
					Filter 1 (nm): 405 Filter 2 (nm): (none) Bicr. factor: 1.00	

Lin limit. low: 0 RBL min (abs): 0.8000
 High: 20000 Max (abs): 2.5000

Calculation. model: Factor RBL stability (days): 1 ** definito dall'operatore
 Factor: 82250.00 Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: NESSUNO
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **CK-NAC FL** - Codes CK F060 / F120 / F245 CH
 Preparation: MIX REAGENTS ACCORDING TO INSERT SHEET

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: CK
 Unit: UI
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	0	0	0	24	204	2000	2000	2000
Female:	0	0	0	24	173	2000	2000	2000
Children:	0	0	0	24	204	2000	2000	2000
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Reaction type: Kinetic
 Direction: Up
 E.P. limit (abs): N/A
 Depl.limit (abs): 2.0000
 First limit (abs): N/A
 Linear factor: N/A
 Fit: 0.970

		Parameter					
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):			0	0	0	170	126
Dil./Rgt. Code:			CK	0			
Lot Number:							
Ratio/Vol (µl):	1/1		300	0	0		
Rinse (µl):			0	0	0		
Sample (µl):			12				

Filter 1 (nm): 340
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 0 RBL min (abs): -0.2000
 High: 2000 Max (abs): 0.8000

Calculation. model: Factor RBL stability (days): 1 ** definito dall'operatore
 Factor: 4127.00 Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: NESSUNO
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **CK-MB FL** - Code MB F060 / F120 CH
 Preparation: MIX REAGENTS ACCORDING TO INSERT SHEET

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: CKMB
 Unit: UI
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	0	0	0	0	24	1200	1200	1200
Female:	0	0	0	0	24	1200	1200	1200
Children:	0	0	0	0	24	1200	1200	1200
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Reaction type: Kinetic
 Direction: Up
 E.P. limit (abs): N/A
 Depl.limit (abs): 2.0000
 First limit (abs): N/A
 Linear factor: N/A
 Fit: 0.970

		Parameter								
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read			
Times (sec):			0	0	0	242	126			
Dil./Rgt. Code:			CKMB	0		<table border="1"> <tr> <td>Filter 1 (nm): 340</td> </tr> <tr> <td>Filter 2 (nm): (none)</td> </tr> <tr> <td>Bicr. factor: 1.00</td> </tr> </table>		Filter 1 (nm): 340	Filter 2 (nm): (none)	Bicr. factor: 1.00
Filter 1 (nm): 340										
Filter 2 (nm): (none)										
Bicr. factor: 1.00										
Lot Number:										
Ratio/Vol (µl):	1/1		300	0	0					
Rinse (µl):			0	0	0					
Sample (µl):			12							

Lin limit. low: 0 RBL min (abs): -0.2000
 High: 1600 Max (abs): 2.0000

Calculation. model: Factor RBL stability (days): 1 ** definito dall'operatore
 Factor: 8254.00 Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: NESSUNO
 Controlli suggeriti: CKMB

Analyzer: ILab 300

Application: **COPPER** - Code CU 0100 CH
 Preparation: MIX REAGENTS ACCORDING TO INSERT SHEET
 Storage: KEEP AT ROOM TEMP (MANDATORY!)
 Stability: 14 DAYS

PROGRAM

Description: COPPER
 Unit: µg/dl
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	5	5	5	70	140	500	500	500
Female:	5	5	5	80	155	500	500	500
Children:	2	2	2	30	150	500	500	500
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Reaction type: End Point
 Direction: Up
 E.P. limit (abs): 0.1000
 Depl.limit (abs): N/A
 First limit (abs): N/A
 Linear factor: N/A
 Fit: N/A

		Parameter					
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):			0	0	0	350	0
Dil./Rgt. Code:			COPPER	0			
Lot Number:							
Ratio/Vol (µl):	1/1		300	0	0		
Rinse (µl):			0	0	0		
Sample (µl):			20				

Filter 1 (nm): 578
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 2 RBL min (abs): -0.2000
 High: 500 Max (abs): 2.0000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 1 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: STANDARD NEL KIT
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **CREATININE** - Code CR 0500 / 1500 CH

Preparation: REAGENT # 1 - USE REAGENT A

REAGENT # 2 - USE REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: UNTIL EXPIRATION DATE ON LABEL

Description: CREATININE

Unit: mg/dl

Decimals: 2

LIS Code: **

Unit Factor: 1.0

Slope: 1.00

Intercept: 0.00

Reference Range

	LOW VALUES				HIGH VALUES			
Male:	0.02	0.02	0.02	0.60	1.20	12.00	12.00	12.00
Female:	0.02	0.02	0.02	0.50	1.00	12.00	12.00	12.00
Children:	0.02	0.02	0.02	0.60	1.20	12.00	12.00	12.00
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No				No	No	No

Parameter

Reaction type: Fixed time

Direction: Up

E.P. limit (abs): N/A

Depl.limit (abs): 2.0000

First limit (abs): 1.0000

Linear factor: 0.98

Fit: N/A

Times (sec):

Dil./Rgt. Code:

Lot Number:

Ratio/Vol (µl):

Rinse (µl):

Sample (µl):

Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
	0	72	0	26	72
	CA1	CA2			
1/5	150	150	0		
	0	0	0		
	20				

Filter 1 (nm): 510
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 0.02

High: 12.00

RBL min (abs): -0.2000

Max (abs): 2.0000

Calculation. model: Standard

Factor: 0.00

Sample blank: N/A

RBL stability (days): 1

Calibration stab. (days): 3

Dynamic controls (min): None

** definito dall'operatore

N/A non applicabile

Calibratore suggerito: AUTOCAL H

Controlli suggeriti: QUANTINORM CHEMA

QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **GAMMA-GT FL** - Codes GT F080 / F245 / F400 / F600 CH
 Preparation: MIX REAGENTS ACCORDING TO INSERT SHEET

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: GAMMA GT

Unit: UI
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	0	0	0	0	50	800	800	800
Female:	0	0	0	0	30	800	800	800
Children:	0	0	0	0	50	800	800	800
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Parameter

Reaction type: Kinetic
 Direction: Up
 E.P. limit (abs): N/A
 Depl.limit (abs): 2.0000
 First limit (abs): N/A
 Linear factor: N/A
 Fit: 0.970

	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read			
Times (sec):		0	0	0	44	144			
Dil./Rgt. Code:		GGT	0		<table border="1"> <tr> <td>Filter 1 (nm): 405</td> </tr> <tr> <td>Filter 2 (nm): (none)</td> </tr> <tr> <td>Bicr. factor: 1.00</td> </tr> </table>		Filter 1 (nm): 405	Filter 2 (nm): (none)	Bicr. factor: 1.00
Filter 1 (nm): 405									
Filter 2 (nm): (none)									
Bicr. factor: 1.00									
Lot Number:									
Ratio/Vol (µl):	1/1	300	0	0					
Rinse (µl):		0	0	0					
Sample (µl):		30							

Lin limit. low: 0 RBL min (abs): -0.2000
 High: 800 Max (abs): 2.0000

Calculation. model: Factor RBL stability (days): 1 ** definito dall'operatore
 Factor: 1280.00 Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: NESSUNO
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **GLUCOSE FL** - Codes GL F400 / 100F / 150F CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: GLUCOSE
 Unit: mg/dl
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	5	5	5	60	110	500	500	500
Female:	5	5	5	60	110	500	500	500
Children:	5	5	5	60	110	500	500	500
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No				No	No	

Reaction type: End Point
 Direction: Up
 E.P. limit (abs): 0.1000
 Depl.limit (abs): N/A
 First limit (abs): N/A
 Linear factor: N/A
 Fit: N/A

		Parameter					
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):			0	0	0	512	0
Dil./Rgt. Code:			GLUC	0			
Lot Number:							
Ratio/Vol (µl):	1/5		360	0	0		
Rinse (µl):			0	0	0		
Sample (µl):			4				

Filter 1 (nm): 510
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 5 RBL min (abs): -0.2000
 High: 550 Max (abs): 0.8000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 10 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: AUTOCAL H
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **GLUCOSE UV FL - Code GL F601 CH**
 Preparation: **MIX REAGENTS ACCORDING TO INSERT SHEET**

Storage: **REFRIGERATE AT 2-8°C**
 Stability: **UNTIL EXPIRATION DATE ON LABEL**

Description: GLUCOSE
 Unit: mg/dl
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	5	5	5	60	110	800	800	800
Female:	5	5	5	60	110	800	800	800
Children:	5	5	5	60	110	800	800	800
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Reaction type: End Point
 Direction: Up
 E.P. limit (abs): 0.1000
 Depl.limit (abs): N/A
 First limit (abs): N/A
 Linear factor: N/A
 Fit: N/A

		Parameter					
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):			0	0	0	242	0
Dil./Rgt. Code:			GLUC	0			
Lot Number:							
Ratio/Vol (µl):	1/5		360	0	0		
Rinse (µl):			0	0	0		
Sample (µl):			4				

Filter 1 (nm): 340
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 5 RBL min (abs): -0.2000
 High: 1000 Max (abs): 0.8000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 10 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: AUTOCAL H
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **GOT/AST FL** - Codes GO F080 / F245 / F400 / F600 CH
 Preparation: MIX REAGENTS ACCORDING TO INSERT SHEET

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: GOT/AST

Unit: UI
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	0	0	0	0	35	400	400	400
Female:	0	0	0	0	31	400	400	400
Children:	0	0	0	0	35	400	400	400
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Parameter

Reaction type: Kinetic
 Direction: Down
 E.P. limit (abs): N/A
 Depl.limit (abs): 2.0000
 First limit (abs): N/A
 Linear factor: N/A
 Fit: 0.960

	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read			
Times (sec):		0	0	0	170	126			
Dil./Rgt. Code:		GOT	0		<table border="1"> <tr> <td>Filter 1 (nm): 340</td> </tr> <tr> <td>Filter 2 (nm): (none)</td> </tr> <tr> <td>Bicr. factor: 1.00</td> </tr> </table>		Filter 1 (nm): 340	Filter 2 (nm): (none)	Bicr. factor: 1.00
Filter 1 (nm): 340									
Filter 2 (nm): (none)									
Bicr. factor: 1.00									
Lot Number:									
Ratio/Vol (µl):	1/1	300	0	0					
Rinse (µl):		0	0	0					
Sample (µl):		30							

Lin limit. low: 0 RBL min (abs): 1.0000
 High: 440 Max (abs): 2.5000

Calculation. model: Factor RBL stability (days): 1 ** definito dall'operatore
 Factor: -1746.00 Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: NESSUNO
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **GPT/ALT FL** - Codes GP F080 / F245 / F400 / F600 CH
 Preparation: MIX REAGENTS ACCORDING TO INSERT SHEET

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: GPT/ALT
 Unit: UI
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	0	0	0	0	45	400	400	400
Female:	0	0	0	0	34	400	400	400
Children:	0	0	0	0	45	400	400	400
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Reaction type: Kinetic
 Direction: Down
 E.P. limit (abs): N/A
 Depl.limit (abs): 2.0000
 First limit (abs): N/A
 Linear factor: N/A
 Fit: 0.960

		Parameter								
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read			
Times (sec):			0	0	0	170	126			
Dil./Rgt. Code:			GPT	0		<table border="1"> <tr> <td>Filter 1 (nm): 340</td> </tr> <tr> <td>Filter 2 (nm): (none)</td> </tr> <tr> <td>Bicr. factor: 1.00</td> </tr> </table>		Filter 1 (nm): 340	Filter 2 (nm): (none)	Bicr. factor: 1.00
Filter 1 (nm): 340										
Filter 2 (nm): (none)										
Bicr. factor: 1.00										
Lot Number:										
Ratio/Vol (µl):	1/1		300	0	0					
Rinse (µl):			0	0	0					
Sample (µl):			30							

Lin limit. low: 0 RBL min (abs): 1.0000
 High: 440 Max (abs): 2.5000

Calculation. model: Factor RBL stability (days): 1 ** definito dall'operatore
 Factor: -1746.00 Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: NESSUNO
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **IRON FZ** - Code FE F245 / F400 CH
 Preparation: REAGENT # 1 - USE REAGENT A
 REAGENT # 2 - USE REAGENT B (READ PREPARATION IN MANUAL INSERT SHEET)

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: IRON
 Unit: µg/dl
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	5	5	5	59	158	500	500	500
Female:	5	5	5	37	145	500	500	500
Children:	5	5	5	59	158	500	500	500
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Reaction type: End Point
 Direction: Up
 E.P. limit (abs): 0.5000
 Depl.limit (abs): N/A
 First limit (abs): N/A
 Linear factor: N/A
 Fit: N/A

		Parameter					
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):			0	180	0	368	0
Dil./Rgt. Code:			IRON1	IRON2			
Lot Number:							
Ratio/Vol (µl):	1/1		300	75	0		
Rinse (µl):			0	0	0		
Sample (µl):			75				

Filter 1 (nm): 578
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 0 RBL min (abs): -0.2000
 High: 1000 Max (abs): 1.0000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 1 N/A non applicabile
 Sample blank: Yes Dynamic controls (min): None

Calibratore suggerito: AUTOCAL H
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **LDH FL** - Codes LD F060 / F120 / F245 CH
 Preparation: MIX REAGENTS ACCORDING TO INSERT SHEET

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: LDH-P
 Unit: UI
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	30	30	30	225	450	4000	4000	4000
Female:	30	30	30	225	450	4000	4000	4000
Children:	30	30	30	225	450	4000	4000	4000
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Reaction type: Kinetic
 Direction: Down
 E.P. limit (abs): N/A
 Depl.limit (abs): 2.0000
 First limit (abs): N/A
 Linear factor: N/A
 Fit: 0.970

		Parameter					
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):			0	0	0	62	108
Dil./Rgt. Code:			LDH	0			
Lot Number:							
Ratio/Vol (µl):	1/1		400	0	0		
Rinse (µl):			0	0	0		
Sample (µl):			4				

Filter 1 (nm): 340
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 30 RBL min (abs): 0.8000
 High: 4000 Max (abs): 2.5000

Calculation. model: Factor RBL stability (days): 1 ** definito dall'operatore
 Factor: -16030.00 Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: NESSUNO
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **LIPASE FL** - Code LP F060 CH
 Preparation: REAGENT # 1 - USE REAGENT A
 REAGENT # 2 - USE REAGENT B

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: LIPASE
 Unit: UI
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	0	0	0	0	63	250	250	250
Female:	0	0	0	0	63	250	250	250
Children:	0	0	0	0	63	250	250	250
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Reaction type: Kinetic
 Direction: Up
 E.P. limit (abs): N/A
 Depl.limit (abs): 2.0000
 First limit (abs): N/A
 Linear factor: N/A
 Fit: 0.970

		Parameter								
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read			
Times (sec):			0	252	0	62	170			
Dil./Rgt. Code:			LIP1	LIP2		<table border="1"> <tr> <td>Filter 1 (nm): 578</td> </tr> <tr> <td>Filter 2 (nm): (none)</td> </tr> <tr> <td>Bicr. factor: 1.00</td> </tr> </table>		Filter 1 (nm): 578	Filter 2 (nm): (none)	Bicr. factor: 1.00
Filter 1 (nm): 578										
Filter 2 (nm): (none)										
Bicr. factor: 1.00										
Lot Number:										
Ratio/Vol (µl):	1/1		300	60	0					
Rinse (µl):			0	0	0					
Sample (µl):			3							

Lin limit. low: 0 RBL min (abs): -0.5000
 High: 250 Max (abs): 2.0000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 99 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: FORNITO CON IL KIT
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **MAGNESIUM** - Code MG 0200 / 0500 CH

Preparation: REAGENT # 1 - USE REAGENT B

REAGENT # 2 - USE REAGENT A

Storage: REFRIGERATE AT 2-8°C

Stability: UNTIL EXPIRATION DATE ON LABEL

Description: MAGNESIUM

Unit: meq/l

Decimals: 2

LIS Code: **

Unit Factor: 1.0

Slope: 1.00

Intercept: 0.00

Reference Range

	LOW VALUES			HIGH VALUES				
Male:	0.00	0.00	0.00	1.30	2.10	8.00	8.00	8.00
Female:	0.00	0.00	0.00	1.30	2.10	8.00	8.00	8.00
Children:	0.00	0.00	0.00	1.30	2.10	8.00	8.00	8.00
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No				No	No	

Parameter

Reaction type: End Point

Direction: Up

E.P. limit (abs): 0.1000

Depl.limit (abs): N/A

First limit (abs): N/A

Linear factor: N/A

Fit: N/A

Times (sec):

Dil./Rgt. Code:

Lot Number:

Ratio/Vol (µl):

Rinse (µl):

Sample (µl):

Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
	0	54	0	98	0
	MG1	MG2			
1/1	300	300	0		
	0	0	0		
	6				

Filter 1 (nm): 510
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 0.00

High: 8.00

RBL min (abs): -0.2000

Max (abs): 2.0000

Calculation. model: Standard

Factor: N/A

Sample blank: Yes

RBL stability (days): 1

Calibration stab. (days): 1

Dynamic controls (min): None

** definito dall'operatore

N/A non applicabile

Calibratore suggerito: AUTOCAL H

Controlli suggeriti: QUANTINORM CHEMA

QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **MAGNESIUM XL** - Code MX 0300 / 0500 CH

Preparation: LIQUID REAGENT READY TO USE

Storage: REFRIGERATE AT 2-8°C

Stability: UNTIL EXPIRATION DATE ON LABEL

Description: MAGNESIUM

Unit: meq/l

Decimals: 2

LIS Code: **

Unit Factor: 1.0

Slope: 1.00

Intercept: 0.00

Reference Range

	LOW VALUES			HIGH VALUES				
Male:	0.00	0.00	0.00	1.30	2.10	8.00	8.00	8.00
Female:	0.00	0.00	0.00	1.30	2.10	8.00	8.00	8.00
Children:	0.00	0.00	0.00	1.30	2.10	8.00	8.00	8.00
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No				No	No	

Parameter

Reaction type: End Point

Direction: Up

E.P. limit (abs): 0.1000

Depl.limit (abs): N/A

First limit (abs): N/A

Linear factor: N/A

Fit: N/A

Times (sec):

Dil./Rgt. Code:

Lot Number:

Ratio/Vol (µl):

Rinse (µl):

Sample (µl):

Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
	0	0	0	98	0
	MG	0			
1/1	360	0	0		
	0	0	0		
	3				

Filter 1 (nm): 546
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 0.00

High: 8.00

RBL min (abs): -0.2000

Max (abs): 2.0000

Calculation. model: Standard

Factor: N/A

Sample blank: Yes

RBL stability (days): 1

Calibration stab. (days): 1

Dynamic controls (min): None

** definito dall'operatore

N/A non applicabile

Calibratore suggerito: AUTOCAL H
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **PHOSPHORUS** - Code PH 0100 / 0500 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: PHOSPHORUS
 Unit: mg/dl
 Decimals: 2
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

Reference Range

	LOW VALUES				HIGH VALUES			
Male:	0.00	0.00	0.00	4.00	6.50	15.00	15.00	15.00
Female:	0.00	0.00	0.00	4.00	6.50	15.00	15.00	15.00
Children:	0.00	0.00	0.00	3.00	4.50	15.00	15.00	15.00
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No				No	No	

Parameter

Reaction type: End Point	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Direction: Up		0	0	0	350	0
E.P. limit (abs): 0.1000		PH	0			
DepL.limit (abs): N/A		Lot Number:				
First limit (abs): N/A	1/5	Ratio/Vol (µl): 400	0	0	<div style="border: 1px solid black; padding: 5px;"> Filter 1 (nm): 340 Filter 2 (nm): (none) Bicr. factor: 1.00 </div>	
Linear factor: N/A		Rinse (µl): 0	0	0		
Fit: N/A		Sample (µl): 4				

Lin limit. low: 0.10 RBL min (abs): -0.2000
 High: 20.00 Max (abs): 2.0000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 7 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: AUTOCAL H
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **PROTEINS (TOTAL)** - Codes TP 0100 / 0500 / 1000 / 1500 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: ROOM TEMPERATURE (2-30°C)
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: PROTEINS
 Unit: g/dl
 Decimals: 2
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

Reference Range

	LOW VALUES			HIGH VALUES				
Male:	1.00	1.00	1.00	6.30	8.40	12.00	12.00	12.00
Female:	1.00	1.00	1.00	6.30	8.40	12.00	12.00	12.00
Children:	1.00	1.00	1.00	6.30	8.40	12.00	12.00	12.00
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No				No	No	

Parameter

Reaction type: End Point	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Direction: Up		0	0	0	548	0
E.P. limit (abs): 0.1000		TP	0			
Depl.limit (abs): N/A		Lot Number:				
First limit (abs): N/A	1/3	Ratio/Vol (µl): 400	0	0	Filter 1 (nm): 546 Filter 2 (nm): (none) Bicr. factor: 1.00	
Linear factor: N/A		Rinse (µl): 0	0	0		
Fit: N/A		Sample (µl): 4				

Lin limit. low: 0.10 RBL min (abs): -0.2000
 High: 12.00 Max (abs): 2.0000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 10 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: AUTOCAL H
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **PROTEINS HS** - Code HS 0100 / 0500 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: ROOM TEMPERATURE (2-30°C)
 Stability: UNTIL EXPIRATION DATE ON LABEL

		Reference Range							
Description: HS PROTEINS		LOW VALUES		HIGH VALUES					
Unit: mg/dl	Male: 2	2	2	28	141	500	500	500	
Decimals: 0	Female: 2	2	2	28	141	500	500	500	
LIS Code: **	Children: 2	2	2	28	141	500	500	500	
Unit Factor: 1.0	Low Alert: Low alert	Very low	Low	Normal values		High	Very high	High alert	
Slope: 1.00	Rerun: No	No				No	No	No	
Intercept: 0.00									

		Parameter						
Reaction type: End Point	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read		
Direction: Up	Times (sec):	0	0	0	548	0		
E.P. limit (abs): 0.1000	Dil./Rgt. Code:	HS	0					
Depl.limit (abs): N/A	Lot Number:				Filter 1 (nm): 620			
First limit (abs): N/A	Ratio/Vol (µl): 1/5	450	0	0	Filter 2 (nm): (none)			
Linear factor: N/A	Rinse (µl):	0	0	0	Bicr. factor: 1.00			
Fit: N/A	Sample (µl):	3						

Lin limit. low: 2 RBL min (abs): -0.2000
 High: 500 Max (abs): 2.0000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 5 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: STANDARD NEL KIT
 Controlli suggeriti: XX

Analyzer: ILab 300

Application: **TRIGLYCERIDES FL** - Codes TR F100 / F400 / 100F / 150F CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: TRIGLYCERIDES

Unit: mg/dl
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	2	2	2	10	190	800	800	800
Female:	2	2	2	10	190	800	800	800
Children:	2	2	2	10	190	800	800	800
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Parameter

Reaction type: End Point
 Direction: Up
 E.P. limit (abs): 0.1000
 Depl.limit (abs): N/A
 First limit (abs): N/A
 Linear factor: N/A
 Fit: N/A

	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):		0	0	0	440	0
Dil./Rgt. Code:		TRIG	0			
Lot Number:						
Ratio/Vol (µl):	1/3	360	0	0		
Rinse (µl):		0	0	0		
Sample (µl):		4				

Filter 1 (nm): 510
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 5 RBL min (abs): -0.2000
 High: 1000 Max (abs): 0.8000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 10 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: AUTOCAL H
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **UREA UV FL** - Codes AZ F080 / F245 / F400 / F600 CH
 Preparation: MIX REAGENTS ACCORDING TO INSERT SHEET

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: UREA
 Unit: mg/dl
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	1	1	1	10	50	300	300	300
Female:	1	1	1	10	50	300	300	300
Children:	1	1	1	10	50	300	300	300
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Reaction type: Fixed time
 Direction: Down
 E.P. limit (abs): N/A
 Depl.limit (abs): 0.3000
 First limit (abs): 1.0000
 Linear factor: 0.95
 Fit: N/A

		Parameter					
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):			0	0	0	26	72
Dil./Rgt. Code:			UREA	0			
Lot Number:							
Ratio/Vol (µl):	1/3		400	0	0		
Rinse (µl):			0	0	0		
Sample (µl):			4				

Filter 1 (nm): 340
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 3 RBL min (abs): 1.0000
 High: 300 Max (abs): 2.0000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: 0.00 Calibration stab. (days): 10 N/A non applicabile
 Sample blank: N/A Dynamic controls (min): None

Calibratore suggerito: AUTOCAL H
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **URIC ACID T FL** - Code AU F402 CH
 Preparation: MIX REAGENTS ACCORDING TO INSERT SHEET

Storage: REFRIGERATE AT 2-8°C
 Stability: UNTIL EXPIRATION DATE ON LABEL

Description: URIC ACID				Reference Range				
Unit: mg/dl	Male: 1.00	1.00	1.00	3.50	7.00	25.00	25.00	25.00
Decimals: 2	Female: 1.00	1.00	1.00	2.40	5.70	25.00	25.00	25.00
LIS Code: **	Children:							
Unit Factor: 1.0	Low Alert: Low alert	Very low	Low	Normal values		High	Very high	High alert
Slope: 1.00	Rerun: No	No				No	No	No
Intercept: 0.00								

Reaction type: End Point	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Direction: Up	Times (sec):	0	0	0	350	0
E.P. limit (abs): 0.1000	Dil./Rgt. Code:	URIC	0			
Depl.limit (abs): N/A	Lot Number:				Filter 1 (nm): 546	
First limit (abs): N/A	Ratio/Vol (µl): 1/5	320	0	0	Filter 2 (nm): (none)	
Linear factor: N/A	Rinse (µl):	0	0	0	Bicr. factor: 1.00	
Fit: N/A	Sample (µl):	10				

Lin limit. low: 0.50 RBL min (abs): -0.2000
 High: 25.00 Max (abs): 0.8000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 10 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: AUTOCAL H
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **URIC ACID AOX FL** - Code AX F100 / F250 / F600 CH

Preparation: REAGENT # 1 - USE REAGENT A

REAGENT # 2 - USE REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: UNTIL EXPIRATION DATE ON LABEL

Description: URIC ACID

Unit: mg/dl

Decimals: 2

LIS Code: **

Unit Factor: 1.0

Slope: 1.00

Intercept: 0.00

Reference Range

	LOW VALUES			HIGH VALUES				
Male:	1.00	1.00	1.00	3.50	7.00	25.00	25.00	25.00
Female:	1.00	1.00	1.00	2.40	5.70	25.00	25.00	25.00
Children:	0	0	0	0	0	0	0	0
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No				No	No	No

Parameter

	Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Reaction type: End Point		0	180	0	368	0
Direction: Up		URIC1	URIC2			
E.P. limit (abs): 0.1000						
DepL.limit (abs): N/A						
First limit (abs): N/A	1/1	300	75	0		
Linear factor: N/A		0	0	0		
Fit: N/A		20				
					Filter 1 (nm): 546 Filter 2 (nm): (none) Bicr. factor: 1.00	

Lin limit. low: 0

High: 35

RBL min (abs): -0.2000

Max (abs): 1.0000

Calculation. model: Standard

Factor: N/A

Sample blank: Yes

RBL stability (days): 1

Calibration stab. (days): 1

Dynamic controls (min): None

** definito dall'operatore

N/A non applicabile

Calibratore suggerito: AUTOCAL H

Controlli suggeriti: QUANTINORM CHEMA

QUANTIPATH CHEMA

Analyzer: ILab 300

Application: **ZINC** - Code ZN 0125 CH
 Preparation: PREPARE MIXED REAGENT AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 14 DAYS

Description: ZINK
 Unit: µg/dl
 Decimals: 0
 LIS Code: **
 Unit Factor: 1.0
 Slope: 1.00
 Intercept: 0.00

		Reference Range						
		LOW VALUES			HIGH VALUES			
Male:	5	5	5	70	150	500	500	500
Female:	5	5	5	70	150	500	500	500
Children:	5	5	5	70	150	500	500	500
Low Alert:	Low alert	Very low	Low	Normal values		High	Very high	High alert
Rerun:	No	No					No	No

Reaction type: End Point
 Direction: Up
 E.P. limit (abs): 0.1000
 Depl.limit (abs): N/A
 First limit (abs): N/A
 Linear factor: N/A
 Fit: N/A

		Parameter					
		Predilut.->	S.+R. 1->	Reag. 2->	Reag. 3->	Incubation ->	Read
Times (sec):			0	0	0	350	0
Dil./Rgt. Code:			ZINC	0			
Lot Number:							
Ratio/Vol (µl):	1/1		400	0	0		
Rinse (µl):			0	0	0		
Sample (µl):			20				

Filter 1 (nm): 578
 Filter 2 (nm): (none)
 Bicr. factor: 1.00

Lin limit. low: 2 RBL min (abs): -0.2000
 High: 1000 Max (abs): 2.0000

Calculation. model: Standard RBL stability (days): 1 ** definito dall'operatore
 Factor: N/A Calibration stab. (days): 1 N/A non applicabile
 Sample blank: No Dynamic controls (min): None

Calibratore suggerito: STANDARD NEL KIT
 Controlli suggeriti: QUANTINORM CHEMA
 QUANTIPATH CHEMA