



ALCYON 300 - FALCOR 300 instructions

rev 431.0.2 - 2013-01-10

Chema
D I A G N O S T I C A

Chema Diagnostica
Via Campania 2/4
I-60030 Monsano | An Italy
T +39 0731 605064 | F +39 0731 605672
www.chema.com | mail@chema.com

ALCYON 300 - FALCOR 300 instructions

rev 431.0.2 - 2013-01-10

Warning!

These programs should be used only as guideline. To check the reagent correct setting on the instrument, use only good quality control sera and work accordingly to with good laboratory practice. Instrumental factor must be checked. Chema Diagnostica will not be responsible for instrumental bad programming.

Analyzer: ALCYON 300 - FALCOR 300

Application: **ALBUMIN** - Code BC 0100 - 0500 - 1500 CH
 Preparation: R1 - LIQUID READY TO USE SINGLE REAGENT
 Storage: ROOM TEMPERATURE OR REFRIGERATE (2-30°C)
 Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	g/dl
NUMBER OF DECIMALS:	1
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	400
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	2.0
LOWER NORMAL VALUE:	3.8
HIGHER NORMAL VALUE:	5.1
DILUTION FACTOR:	2
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.3000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	600
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.050
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	2:00
INCUBATION:	0:00
TOTAL:	1:00

Analyzer: ALCYON 300 - FALCOR 300Application: **ALKALINE PHOSPHATASE FL DGKC** - Codes AL F080 - F245 - F400 - F600 CH

Preparation: MIX 4 PARTS OF REAGENT A + 1 PART OF REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: 30 DAYS AT 2-8°C

PROGRAM

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	4.0
LOWER NORMAL VALUE:	64
HIGHER NORMAL VALUE:	306
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.4000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	405
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.3500
C COMPENSATION FACTOR:	0.050
CALIBRATION FACTOR:	5514
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	2:00
INCUBATION:	0:00
TOTAL:	2:00

Analyzer: ALCYON 300 - FALCOR 300Application: **ALKALINE PHOSPHATASE FL IFCC** - Codes AF F080 - F245 - F400 CH

Preparation: MIX 4 PARTS OF REAGENT A + 1 PART OF REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: 30 DAYS AT 2-8°C

PROGRAM

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	4.0
LOWER NORMAL VALUE:	35
HIGHER NORMAL VALUE:	104
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.4000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	405
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.3500
C COMPENSATION FACTOR:	0.050
CALIBRATION FACTOR:	5514
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	2:00
INCUBATION:	0:00
TOTAL:	2:00

Analyzer: ALCYON 300 - FALCOR 300

Application: **AMYLASE FL** - Codes AM F060 - F120 - F245
 Preparation: R1 - LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	5.0
LOWER NORMAL VALUE:	0
HIGHER NORMAL VALUE:	96
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.2000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	405
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.3500
C COMPENSATION FACTOR:	0.050
CALIBRATION FACTOR:	6350
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	1:00
INCUBATION:	0:00
TOTAL:	1:00

Analyzer: ALCYON 300 - FALCOR 300

Application: **AMYLASE EPS FL** - Codes EA F080 - F245 CH
Preparation: MIX 4 PARTS OF REAGENT A + 1 PART OF REAGENT B
Storage: REFRIGERATE AT 2-8°C
Stability: 30 DAYS AT 2-8°C.

PROGRAM

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	4.0
LOWER NORMAL VALUE:	10
HIGHER NORMAL VALUE:	220
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.2000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	405
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.3500
C COMPENSATION FACTOR:	0.050
CALIBRATION FACTOR:	23200
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	3:00
INCUBATION:	0:00
TOTAL:	3:00

Analyzer: ALCYON 300 - FALCOR 300

Application: **PANCREATIC ISOAMYLASE FL** - Codes PA F080 CH

Preparation: R1 - INSTALL LIQUID READY TO USE REAGENT A
 R2 - INSTALL LIQUID READY TO USE REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	50
SAMPLE VOLUME:	4.0
LOWER NORMAL VALUE:	17
HIGHER NORMAL VALUE:	115
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.2000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	405
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.3500
C COMPENSATION FACTOR:	0.050
CALIBRATION FACTOR:	30280
SLOPE:	1.00
INTERCEPT:	0.00
DELAY FOR RGT. 2 DISPENSE:	3:00
REACTION DELAY:	2:00

Analyzer: ALCYON 300 - FALCOR 300

Application: **BICARBONATE FL** - Codes BR F060 / F245 / F400 CH
 Preparation: R1 - LIQUID READY TO USE SINGLE REAGENT
 Storage: ROOM TEMPERATURE OR REFRIGERATE (2-30°C)
 Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	mmol/l
NUMBER OF DECIMALS:	1
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	300
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	22
HIGHER NORMAL VALUE:	29
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.3000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	405
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.1000
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	5:12

Analyzer: ALCYON 300 - FALCOR 300Application: **BILIRUBIN DIRECT FL** - Code DD F125 - F500 CHPreparation: R1 - INSTALL LIQUID READY TO USE REAGENT A
R2 - INSTALL LIQUID READY TO USE REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	1
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	10.0
LOWER NORMAL VALUE:	0.0
HIGHER NORMAL VALUE:	0.2
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.0800
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	550
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	YES
T TOLERANCE FACTOR:	0.050
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	2:12

Analyzer: ALCYON 300 - FALCOR 300

Application: **BILIRUBIN TOTAL FL** - Code DT F125 - F500 CH
Preparation: R1 - INSTALL LIQUID READY TO USE REAGENT A
R2 - INSTALL LIQUID READY TO USE REAGENT B
Storage: REFRIGERATE AT 2-8°C
Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	1
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	50
SAMPLE VOLUME:	10.0
LOWER NORMAL VALUE:	0.1
HIGHER NORMAL VALUE:	1.2
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.0800
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	500
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
T TOLERANCE FACTOR:	0.050
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
DELAY FOR RGT.2 DISPENSE:	0:36
REACTION DELAY:	5:12

Analyzer: ALCYON 300 - FALCOR 300
Application: **CALCIUM** - Code CA 0300 - 0505 CH
Preparation: MIX 1 PART OF REAGENT A + 1 PART OF REAGENT B
Storage: REFRIGERATE AT 2-8°C
Stability: 14 DAYS

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	1
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	250
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	4.0
LOWER NORMAL VALUE:	8.8
HIGHER NORMAL VALUE:	10.8
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.2500
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	550
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.050
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	2:12

Analyzer: ALCYON 300 - FALCOR 300
 Application: **CHOLESTEROL FL** - Code CT F100 - F400 - 100F CH
 Preparation: R1 - LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	300
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	140
HIGHER NORMAL VALUE:	200
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.1000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	500
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.050
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	4:48

Analyzer: ALCYON 300 - FALCOR 300Application: **HDL DIRECT FL-** Code HD F080 - F245 - F400 CHPreparation: R1 - INSTALL LIQUID READY TO USE REAGENT A
R2 - INSTALL LIQUID READY TO USE REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: 60 DAYS OPEN BOTTLES- CALIBRATOR: AS INDICATED IN INSERT SHEET

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	210
SECOND REAGENT VOL.:	70
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	35
HIGHER NORMAL VALUE:	55
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.3000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	600
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.050
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
DELAY FOR RGT. 2 DISPENSE:	3:00
REACTION DELAY:	3:24

Analyzer: ALCYON 300 - FALCOR 300

Application: **LDL DIRECT FL** - Code DL F080 CH
 Preparation: R1 - USE LIQUID READY TO USE REAGENT A
 R2 - USE LIQUID READY TO USE REAGENT B
 Storage: REFRIGERATE AT 2-8°C
 Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	210
SECOND REAGENT VOL.:	70
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	76
HIGHER NORMAL VALUE:	218
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.3000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	600
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.050
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
DELAY FOR RGT. 2 DISPENSE:	3:00
REACTION DELAY:	3:24

Analyzer: ALCYON 300 - FALCOR 300Application: **CHOLINESTERASE FL (DGKC)** - Codes CH F096 / F245 CH

Preparation: MIX 5 PARTS OF REAGENT A + 1 PART OF REAGENT B

Storage: REFRIGERATE AT 2-8°C

Stability: 14 DAYS AT 2-8°C

PROGRAM

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	240
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	4.0
LOWER NORMAL VALUE:	5600
HIGHER NORMAL VALUE:	11200
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.4000
HIGH REAGENT BLANK:	3.0000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	405
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.3500
C COMPENSATION FACTOR:	0.050
CALIBRATION FACTOR:	131600
SLOPE:	2.00
INTERCEPT:	0.00
REACTION DELAY:	1:36

Analyzer: ALCYON 300 - FALCOR 300
Application: **CK-MB FL** - Code MB F060 - F120 CH
Preparation: MIX 4 PARTS OF REAGENT A + 1 PART OF REAGENT B
Storage: REFRIGERATE AT 2-8°C
Stability: 10 DAYS

PROGRAM

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	10.0
LOWER NORMAL VALUE:	0
HIGHER NORMAL VALUE:	24
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.5000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	340
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.2000
C COMPENSATION FACTOR:	0.1000
CALIBRATION FACTOR:	16500
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	3:00

Analyzer: ALCYON 300 - FALCOR 300
Application: **CK-NAC FL** - Codes CK F060 - F120 - F245 CH
Preparation: **MIX 4 PARTS OF REAGENT A + 1 PART OF REAGENT B**
Storage: **REFRIGERATE AT 2-8°C**
Stability: **14 DAYS**

PROGRAM

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	5.0
LOWER NORMAL VALUE:	0
HIGHER NORMAL VALUE:	190
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.5000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	340
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.2000
C COMPENSATION FACTOR:	0.1000
CALIBRATION FACTOR:	13000
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	1:00

Analyzer: ALCYON 300 - FALCOR 300
Application: **COPPER** - Code CU 0100 CH
Preparation: MIX 1 PART OF REAGENT A + 1 PART OF REAGENT B
Storage: REFRIGERATE AT 2-8°C
Stability: 14 DAYS (WARNING: REAGENT MAY CONDENSE)

PROGRAM

TYPE:	END PT.
UNIT:	microg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	15.0
LOWER NORMAL VALUE:	80
HIGHER NORMAL VALUE:	140
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.3000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	600
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.1000
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	4:24

Analyzer: ALCYON 300 - FALCOR 300
 Application: **CREATININE** - Code CR 0500 - 1000 CH
 Preparation: MIX 1 PART OF REAGENT A + 1 PART OF REAGENT B
 Storage: ROOM TEMPERATURE OR REFRIGERATE (2-30°C)
 Stability: 7 DAYS

PROGRAM

TYPE:	KIN 1
UNIT:	mg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	10.0
LOWER NORMAL VALUE:	0.6
HIGHER NORMAL VALUE:	1.4
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.5000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	500
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
DEPLETION LIMIT:	0.4000
TOLERANCE FACTOR:	0.2000
C COMPENSATION FACTOR:	1.0000
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	0:36
INCUBATION	1:00
TOTAL	1:36

Analyzer: **ALCYON 300 - FALCOR 300**
Application: **CREATININE E FL - Code CE F125 - F375 - F600 CH**
Preparation: R1 - INSTALL LIQUID READY TO USE REAGENT A
 R2 - INSTALL LIQUID READY TO USE REAGENT B
Storage: REFRIGERATE AT 2-8°C
Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	KIN 1
UNIT:	mg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	50
SAMPLE VOLUME:	10.0
LOWER NORMAL VALUE:	0.6
HIGHER NORMAL VALUE:	1.2
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.5000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	500
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
DEPLETION LIMIT:	0.4000
TOLERANCE FACTOR:	0.2000
C COMPENSATION FACTOR:	1.2000
SLOPE:	1.00
INTERCEPT:	0.00
DELAY FOR RGT.2 DISPENSE:	3:24
REACTION DELAY:	1:24
INCUBATION	1:24
TOTAL	6:12

Analyzer: ALCYON 300 - FALCOR 300

Application: **GAMMA-GT FL** - Codes GT F080 - F245 - F400 - F600 CH
Preparation: MIX 4 PARTS OF REAGENT A + 1 PART OF REAGENT B
Storage: REFRIGERATE AT 2-8°C
Stability: 60 DAYS

PROGRAM

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	20.0
LOWER NORMAL VALUE:	10
HIGHER NORMAL VALUE:	47
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.5000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	405
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.3500
C COMPENSATION FACTOR:	0.050
CALIBRATION FACTOR:	2316
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	2:00

Analyzer: ALCYON 300 - FALCOR 300

Application: **GLUCOSE FL** - Code GL F400 - 100F CH
 Preparation: R1 - LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	300
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	60
HIGHER NORMAL VALUE:	110
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.1000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	500
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.1000
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	5:12

Analyzer: ALCYON 300 - FALCOR 300
Application: **GLUCOSE UV FL** - Code GL F601 CH
Preparation: MIX 4 PARTS OF REAGENT A + 1 PART OF REAGENT B
Storage: REFRIGERATE AT 2-8°C
Stability: 60 DAYS

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	300
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	60
HIGHER NORMAL VALUE:	110
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.3000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	340
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.1000
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	5:12

Analyzer: ALCYON 300 - FALCOR 300

Application: **GOT/AST FL** - Codes GO F080 - F245 - F400 - F600 CH
 Preparation: MIX 4 PARTS OF REAGENT A + 1 PART OF REAGENT B
 Storage: REFRIGERATE AT 2-8°C
 Stability: 30 DAYS

PROGRAM

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	20.0
LOWER NORMAL VALUE:	0
HIGHER NORMAL VALUE:	40
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.4000
HIGH REAGENT BLANK:	3.0000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	340
HIGH WAVELENGTH:	380
LOW WAVELENGTH:	---
POLYCHROMATIC FACTOR:	0.1500
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.3500
C COMPENSATION FACTOR:	0.3000
CALIBRATION FACTOR:	3700
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	3:00

Analyzer: ALCYON 300 - FALCOR 300Application: **GPT/ALT FL** - Codes GP F080 - F245 - F400 - F600 CHPreparation: **MIX 4 PARTS OF REAGENT A + 1 PART OF REAGENT B**Storage: **REFRIGERATE AT 2-8°C**Stability: **30 DAYS****PROGRAM**

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	20.0
LOWER NORMAL VALUE:	0
HIGHER NORMAL VALUE:	40
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.4000
HIGH REAGENT BLANK:	3.0000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	340
HIGH WAVELENGTH:	380
LOW WAVELENGTH:	---
POLYCHROMATIC FACTOR:	0.1500
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.3500
C COMPENSATION FACTOR:	0.3000
CALIBRATION FACTOR:	3700
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	3:00

Analyzer: ALCYON 300 - FALCOR 300Application: **IRON FZ** - Codes FE F245 / F400 CHPreparation: R1 - INSTALL LIQUID READY TO USE REAGENT A
R2 - INSTALL REAGENT B (READ PREPARATION ON INSTRUCTIONS)

Storage: REFRIGERATE AT 2-8°C

Stability: AS INDICATED IN THE INSTRUCTIONS FOR USE

PROGRAM

TYPE:	END PT.
UNIT:	microg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	50
SAMPLE VOLUME:	25.0
LOWER NORMAL VALUE:	59
HIGHER NORMAL VALUE:	158
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.3000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	550
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
T TOLERANCE FACTOR:	0.050
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
DELAY FOR RGT.2 DISPENSE:	3:00
REACTION DELAY:	6:24

Analyzer: ALCYON 300 - FALCOR 300

Application: **LDH FL** - Codes LD F060 - F120 - F245
Preparation: MIX 4 PARTS OF REAGENT A + 1 PART OF REAGENT B
Storage: REFRIGERATE AT 2-8°C
Stability: 30 DAYS

PROGRAM

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	300
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	225
HIGHER NORMAL VALUE:	450
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.4000
HIGH REAGENT BLANK:	3.0000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	340
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.2000
C COMPENSATION FACTOR:	0.1000
CALIBRATION FACTOR:	59400
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	1:00

Analyzer: ALCYON 300 - FALCOR 300

Application: **LIPASE FL** - Code LP F060 CH
Preparation: R1 - INSTALL LIQUID READY TO USE REAGENT A
R2 - INSTALL LIQUID READY TO USE REAGENT B
Storage: REFRIGERATE AT 2-8°C
Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	KIN 0
UNIT:	UI/L
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	240
SECOND REAGENT VOL.:	60
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	0
HIGHER NORMAL VALUE:	63
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.5000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	600
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
HIGH LIMIT:	0.15
DEPLETION LIMIT:	0.3500
C COMPENSATION FACTOR:	0.050
SLOPE:	1.00
INTERCEPT:	0.00
DELAY FOR RGT. 2 DISPENSE:	1:24
REACTION DELAY:	3:00
TOTAL:	4:24

Analyzer: **ALCYON 300 - FALCOR 300**
Application: **MAGNESIUM** - Codes MG 0200 / 0500 CH
Preparation: MIX 1 PART OF REAGENT A + 1 PART OF REAGENT B
Storage: ROOM TEMPERATURE OR REFRIGERATE (2-30°C)
Stability: 90 DAYS

PROGRAM

TYPE:	END PT.
UNIT:	mEq/l
NUMBER OF DECIMALS:	1
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	300
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	1.3
HIGHER NORMAL VALUE:	2.1
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.3000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	500
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.1000
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	2:12

Analyzer: ALCYON 300 - FALCOR 300

Application: **MAGNESIUM XL** - Codes MX 0300 / 0500 CH
 Preparation: R1 - LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE (2-8°C)
 Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	mEq/l
NUMBER OF DECIMALS:	1
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	300
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	1.3
HIGHER NORMAL VALUE:	2.1
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.3000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	550
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.1000
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	2:12

Analyzer: ALCYON 300 - FALCOR 300

Application: **PHOSPHORUS** - Code PH 0100 - 0500 CH
Preparation: R1 - LIQUID READY TO USE SINGLE REAGENT
Storage: REFRIGERATE AT 2-8°C
Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	1
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	300
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	3.0
HIGHER NORMAL VALUE:	4.5
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.3000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	340
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.050
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	4:48

Analyzer: ALCYON 300 - FALCOR 300

Application: **PROTEINS HS** - Code HS 0100 - 0500 CH
 Preparation: R1 - LIQUID READY TO USE SINGLE REAGENT
 Storage: ROOM TEMPERATURE OR REFRIGERATE (2-30°C)
 Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	300
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	28
HIGHER NORMAL VALUE:	141
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.2000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	600
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.1000
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	4:24

Analyzer: ALCYON 300 - FALCOR 300

Application: **TOTAL PROTEINS** - Codes TP 0100 - 0500 - 1000 CH
 Preparation: R1 - LIQUID READY TO USE SINGLE REAGENT
 Storage: ROOM TEMPERATURE OR REFRIGERATE (2-30°C)
 Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	g/dl
NUMBER OF DECIMALS:	1
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	300
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	6.3
HIGHER NORMAL VALUE:	8.4
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.2000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	550
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.1000
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	4:48

Analyzer: ALCYON 300 - FALCOR 300
 Application: **TRIGLYCERIDES FL** - Code TR F100 - F400 - 100F CH
 Preparation: R1 - LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	300
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	10
HIGHER NORMAL VALUE:	190
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.1000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	500
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.050
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	4:48

Analyzer: ALCYON 300 - FALCOR 300

Application: **UREA UV FL** - Code AZ F080 - F245 - F400 - F600 CH
Preparation: MIX 4 PARTS OF REAGENT A + 1 PART OF REAGENT B
Storage: REFRIGERATE AT 2-8°C
Stability: 30 DAYS

PROGRAM

TYPE:	KIN 1
UNIT:	mg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	300
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	3.0
LOWER NORMAL VALUE:	10
HIGHER NORMAL VALUE:	50
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.4000
HIGH REAGENT BLANK:	3.0000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	340
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
DEPLETION LIMIT:	0.4000
TOLERANCE FACTOR:	0.2000
C COMPENSATION FACTOR:	1.0000
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	0:36
INCUBATION	1:00
TOTAL	1:36

Analyzer: ALCYON 300 - FALCOR 300Application: **URIC ACID T FL** - Code AU F250 - F402 CH

Preparation: R1 - MIX REAGENT A + REAGENT B AS INDICATED IN INSERT SHEET

Storage: REFRIGERATE AT 2-8°C

Stability: 90 DAYS

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	1
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	5.0
LOWER NORMAL VALUE:	3.5
HIGHER NORMAL VALUE:	7.0
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.1000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	550
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.050
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	4:48

Analyzer: **ALCYON 300 - FALCOR 300**
Application: **URIC ACID AOX FL - Code AX F125 - F245 - F600 CH**
Preparation: R1 - INSTALL LIQUID READY TO USE REAGENT A
 R2 - INSTALL LIQUID READY TO USE REAGENT B
Storage: REFRIGERATE AT 2-8°C
Stability: AS INDICATED IN THE LABEL

PROGRAM

TYPE:	END PT.
UNIT:	mg/dl
NUMBER OF DECIMALS:	1
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	50
SAMPLE VOLUME:	15.0
LOWER NORMAL VALUE:	3.5
HIGHER NORMAL VALUE:	7.0
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.1000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	550
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.050
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
DELAY FOR RGT.2 DISPENSE:	3:00
REACTION DELAY:	6:24

Analyzer: ALCYON 300 - FALCOR 300
 Application: **ZINC** - Code ZN 0125 CH
 Preparation: MIX 4 PARTS OF REAGENT A + 1 PART OF REAGENT B
 Storage: REFRIGERATE AT 2-8°C
 Stability: 2 WEEKS

PROGRAM

TYPE:	END PT.
UNIT:	microg/dl
NUMBER OF DECIMALS:	0
DILUENT:	97 H2O D
FIRST REAGENT VOL.:	200
SECOND REAGENT VOL.:	0
SAMPLE VOLUME:	10.0
LOWER NORMAL VALUE:	70
HIGHER NORMAL VALUE:	115
DILUTION FACTOR:	4
SYSTEMATIC DILUTION:	1
DILUTION FOR URINE:	1
LOW REAGENT BLANK:	0.0000
HIGH REAGENT BLANK:	0.3000
MAXIMUM DEVIATION:	0.0500
PRIMARY WAVELENGTH:	600
HIGH WAVELENGTH:	---
LOW WAVELENGTH:	---
SYSTEMATIC SERUM BLANK:	NO
T TOLERANCE FACTOR:	0.1000
ANTIGEN EXCESS TEST (%):	NO
SLOPE:	1.00
INTERCEPT:	0.00
REACTION DELAY:	4:24