

# AMYLASE FL

AM F060 CH	6 x 10 ml
AM F120 CH	12 x 10 ml
AM F245 CH	12 x 20 ml

## INTENDED USE

Reagent for quantitative in vitro determination of amylase in biological fluids.

## SUMMARY OF TEST

Assays of amylase activity in serum and urine are largely of use in the diagnosis of diseases of the pancreas and in the investigation of pancreatic function.

## PRINCIPLE OF THE METHOD

The enzyme  $\alpha$ -amylase (EC 3.2.1.1, 1,4  $\alpha$ -D-glucose glucanohydrolase) hydrolyzes the 2-chloro-4-nitrophenyl- $\alpha$ -D-maltotriose (CNP3) to release 2-chloro-4-nitrophenol and form 2-chloro-4-nitrophenyl- $\alpha$ -D-maltoside (CNP2), maltotriose (G3) and glucose (G). The rate of formation of the 2-chloro-4-nitrophenol can be detected spectrophotometrically at 405 nm to give a measurement of  $\alpha$ -amylase activity in the sample.

## KIT COMPONENTS

### For in vitro diagnostic use only.

The components of the kit are stable until expiration date on the label.

Keep away from direct light sources.

### DO NOT PIPETTE BY MOUTH!

**AMY R1**      **F060: 6 x 10 ml (liquid) blue cap**  
**F120: 12 x 10 ml (liquid) blue cap**  
**F245: 12 x 20 ml (liquid) blue cap**

Composition: CNP-G3 2.3 mM, NaCl 350 mM, calcium acetate 6 mM, potassium thiocyanate 600 mM, Good's buffer pH 6.0 100 mM, stabilizers and non-reactive components.

Store all components at 2-8°C.

## MATERIALS REQUIRED BUT NOT SUPPLIED

Current laboratory instrumentation. Spectrophotometer UV/VIS with thermostatic cuvette holder. Automatic micropipettes. Glass or high quality polystyrene cuvettes. Saline solution.

## REAGENT PREPARATION

Use reagent ready to use.

Stability: up to expiration date on labels at 2-8°C.

Stability since first opening of vials: preferably within 60 days at 2-8°C.

## PRECAUTIONS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

## SPECIMEN

Serum, plasma (heparinate only).

Amylase is stable in serum and plasma sample up to 2 months at 2-8°C.

## TEST PROCEDURE

Wavelength:	405 nm
Ligthpath:	1 cm
Temperature:	37°C
dispense in cuvette working reagent:	1 ml
preincubate at 37°C for 5 minutes.	
add sample:	25 $\mu$ l
Mix, execute a first reading of absorbance after 1 minute, incubating at 37°C. Perform other 3 readings at 60 seconds intervals. Calculate the $\Delta A/min$ .	

## RESULTS CALCULATION

Perform calculation in units per litre, multiplying the  $\Delta A/min$  by the factor as it is indicated.

Calculation in U/l:  $\Delta A/min \times 3178$

Calculation in  $\mu$ kat/l:  $U/l \times 0.0167 = \mu$ kat/l

## EXPECTED VALUES

Serum - plasma: < 96 U/l ( $< 1.60 \mu$ kat/l)  
Urine: < 480 U/l ( $< 8.00 \mu$ kat/l)

Each laboratory should establish appropriate reference intervals related to its population.

## QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

### QUANTINORM CHEMA

with normal or close to normal control values

### QUANTIPATH CHEMA

with pathological control values.

If required, a multiparametric, human based calibrator is available:

### AUTOCAL H

Please contact Customer Care for further informations.

## TEST PERFORMANCE

### Linearity

the method is linear up to 2000 U/l.

If a  $\Delta A/min$  of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

### Sensitivity/limit of detection (LOD)

the limit of detection is 0.91 U/l.

### Interferences

no interference was observed by the presence of:

hemoglobin  $\leq 500$  mg/dl  
bilirubin  $\leq 50$  mg/dl  
lipids  $\leq 1200$  mg/dl

### Precision

intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	67.89	0.97	1.42
sample 2	171.67	2.61	1.52

inter-assay (n=21)	mean (U/l)	SD (U/l)	CV%
sample 1	67.81	1.93	2.85
sample 2	175.16	4.92	2.81

### Methods comparison

a comparison between Chema and a commercially available product gave the following results:

Amylase Chema = x  
Amylase competitor = y  
n = 181

$y = 1.071x - 0.54$  U/l     $r^2 = 0.997$

## WASTE DISPOSAL

This product is made to be used in professional laboratories.

P501: Dispose of contents according to national/international regulations.








## REFERENCES

Ranson, JHC. Curr Prob Surg 1979; 16:1.  
Salt WB II, Schenker S. Medicine 1976; 55:269.  
Stefanini P, Ermini M, Carboni M. J Am Surg 1965; 110:866.  
Henry RJ, Chiamori N. Clin Chem 1960; 6:434.  
Kaufman RA, Tielz NW. Clin Chem 1980; 26:846.  
Blair HE. U.S. Patent No. 4.649,108.  
Chavez RG et al. U.S. Patent 4,963,479.  
Demetriou J et al. Clinical Chemistry 1974; Principles and Techniques, 2nd Ed, Harper & Row.  
Young OS, Pestaner LC, Gibberman V. Clin Chem 1975; 21:10.

## MANUFACTURER

Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

## SYMBOLS

	in vitro diagnostic medical device
	batch code
	catalogue number
	temperature limit
	use by date
	caution
	consult instructions for use