# ACIDO URICO T FL

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## DESTINAZIONE D'USO

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## BIBLIOGRAFIA

1. N. Rifai, A.R. Horvath et al. Tiez Textbook of Clinical Chemistry and Molecular Diagnostics, sixth edition 2018.
2. A. Kunst, B. Draeger B. et al. Bergmeyer: Methods of Enzymatic Analysis, third edition 1984.
3. M. Ciaccio, G. Lippi. Biochimica Clinica e Medicina di laboratorio, III edizione 2020, EdiSES Università S.r.l.
4. D.A. Armbruster and T. Pry. Limit of Blank, Limit of Detection and Limit of Quantitation. *Clin Biochem Rev.* 2008; 29(suppl 1): 49-52.
5. M.R. Glick, K.W. Ryder et al. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin. Chem.* 1986; 32: 470-475.
6. E. Theodorsson, B. Magnusson et al. Bias in Clinical Chemistry. *Bioanalysis* 2014; 6(21): 2855-2875.
7. M. Vidali, M. Tronchin et al. Protocollo per la comparazione di due metodi analitici di laboratorio. *Biochimica Clinica* 2016; 40(2): 129-142.
8. CLSI EP0-A:2003 Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.
9. H. Passing and W. Bablok. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. *J. Clin. Chem. Biochem.* 1983; 21: 709-720.
10. L. Bilić-Zulle. Comparison of methods: Passing and Bablok regression. *Biochimica Medica* 2011; 21(1): 49-52.
11. J.R. Genzen, J.J. Hunsaker et al. N-acetylcysteine interference of Trinder-based assays. *Clin Biochem.* 2016; 49(1-2):100-104
12. O. Wiewiorka, Z. Čermáková et al. Drug interference in Trinder reaction. *Euromedlab.* 2017; ISSN 1437-4431
13. D. Barham, P. Trinder. An Improved Colour Reagent for the Determination of Blood Glucose by the Oxidase System. *Analyst.* 1972; 97: 142-145.
14. P. Fossati, L. Prencipe et al. Use of 3, 5-Dichloro-2-hydroxybenzenesulfonic Acid/4-Ami nophenazone Chromogenic System in Direct Enzymic Assay of Uric Acid in Serum and Urine. *Clin. Chem.* 1980; 26(2): 227-231.
15. M. Jelkic-Stankov, P. Djurdjevic et al. Determination of uric acid in human serum by an enzymatic method using N-methyl-N-(4-aminophenyl)-3-methoxyaniline reagent. *J. Serb. Chem. Soc.* 2003; 68 (8-9): 691-698.
16. A. J. Bakker, M. Mücke. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *ClinChemLabMed* 2007;45(9):1240-1243.
17. Bakul I. Dalal, MD, FRCPet et al., Factitious Biochemical Measurements Resulting From Hematologic Conditions. *Am J Clin Pathol* 2009;131:195-204
18. Khandpur, Raghbir Singh. Clinical Chemistry Analyser, Random Access. *Compendium of Biomedical Instrumentation* 2020; 457-460.
19. Data on file

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Chema Diagnostica Srl