APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	PH
APP. CODE:	344
WAVELENGTH (Sec/Pri):	700 - 340
ASSAY:	1-POINT <i>TIME: 10</i> <i>POINT: 31</i>
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5
	R1 VOLUME: 300 R2 VOLUME: 0 R3 VOLUME: 0 R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.00
EXPECTED VALUE:	2.0 - 7.0
UNIT:	mg/dl
INSTR. FACTOR (y=ax+b):	a=1 b= 0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 874)	
TEST NAME:	PH
SAMPLE:	Volume 3 µl Dilution 0 µl
REAGENTS:	R1 Volume 300 µl Dilution 0 µl R2 Volume 0 µl Dilution 0 µl
WAVELENGTH:	Pri. 340 Sec. 700
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 Last 17
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 0.4 H 20
CORRELATION FACTOR:	A 1 B 0
UNIT:	mg/dl
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

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IUS-7.5 AUTO

SPECIMEN	
Serum is the preferred specimen. Although heparinized plasma is acceptable, levels of inorganic phosphate are about 0.2 to 0.3 mg/dl lower than in serum. Anticoagulants such as citrate, oxalate, and EDTA interfere with formation of the phosphomolybdate complex and should not be used. Inorganic phosphate in whole blood specimens may either decrease or increase with time, depending on the type of specimen, temperature, and duration of storage. Levels in plasma or serum are increased by prolonged storage with cells at room temperature or 37 °C; it is important to promptly separate serum or plasma from erythrocytes. Hemolyzed specimens are unacceptable because erythrocytes contain high concentrations of organic phosphate esters, which can be hydrolyzed to inorganic phosphate during storage. Inorganic phosphate increases by 4 to 5 mg/dl per day in hemolyzed specimens stored at 4°C. Glucose phosphate, creatine phosphate, and other organic phosphates may also be hydrolyzed by assay conditions, resulting in overestimation of inorganic phosphate levels. Phosphate is considered to be stable in serum that has been separated from the clot for days at 4°C and months when frozen. Urine samples should be collected in 6 mol/l HCl, 20-30 ml for a 24 hours specimen, to avoid precipitation of phosphate complexes. Dilute urine samples 1:20 with purified water before assay.	
EXPECTED VALUES	
serum/plasma (adults):	2.5 - 4.5 mg/dl (0.81 - 1.45 mmol/l)
serum/plasma (children):	4.0 - 7.0 mg/dl (1.29 - 2.26 mmol/l)
urine (nonrestricted diet):	0.4 - 1.3 g/24h (12.9 - 42.2 mmol/24h)
Each laboratory should establish appropriate reference intervals related to its population.	
QUALITY CONTROL AND CALIBRATION	
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available: QUANTINORM CHEMA - MULTINORM CHEMA with normal or close to normal control values QUANTIPATH CHEMA - MULTIPATH CHEMA with pathological control values. If required, a multiparametric, human based calibrator is available: AUTOCAL H	
Please contact Customer Care for further information.	
TEST PERFORMANCE	
Linearity the method is linear up to 20 mg/dl. If the limit value is exceeded, it is suggested to dilute sample 1+9 with distilled water and to repeat the test, multiplying the result by 10.	
Sensitivity/limit of detection (LOD) the limit of detection is 0.4 mg/dl.	
Interferences no interference was observed by the presence of: bilirubin ≤ 25 mg/dl hemoglobin ≤ 100 mg/dl Hemolysis interferes. Both positive and negative interferences with lipemic samples has been observed.	
Precision	intra-assay (n=10) mean (mg/dl) SD (mg/dl) CV% sample 1 3.50 0.05 1.50 sample 2 5.87 0.11 1.90
	inter-assay (n=20) mean (mg/dl) SD (mg/dl) CV% sample 1 3.41 0.08 2.40 sample 2 5.84 0.11 1.90
Methods comparison a comparison between Chema and a commercially available product gave the following results:	
	phosphorus Chema = x phosphorus competitor = y n = 102 y = 1.005x - 0.109 mg/dl r ² = 0.975
WASTE DISPOSAL	
This product is made to be used in professional laboratories. P501: Dispose of contents according to national/international regulations.	

IUS-7.5 AUTO

4/4

1/4