**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
HITACHI 911/912**

TEST:	AMY		
APP. CODE:	334		
WAVELENGTH (Sec/Pri):	700 - 415		
ASSAY:	RATE-A	TIME:	10
		POINT:	7 - 14
SAMPLE VOL:	NORMAL: 7		
	DECREASE: 5		
	INCREASE: 10		
	R1 VOLUME: 280		
	R2 VOLUME: 0		
	R3 VOLUME: 0		
	R4 VOLUME: 0		
ABS LIMIT:	32000 - INC		
PROZONE LIMIT:	0 - UPPER		
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)		
SD LIMIT:	0.250		
DUPLICATE LIMIT:	3%		
ST. 1 CONC:	0.00		
EXPECTED VALUE:	0 - 96		
UNIT:	U/l		
INSTR. FACTOR (y=ax+b):	a=1	b=0	

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
OLYMPUS AU 400/480/600/640/680/2700 (Test code 859)**

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INTENDED USE

Reagent for quantitative in vitro determination of amylase in biological fluids.

PRINCIPLE OF THE METHOD

The enzyme α-amylase hydrolyzes the 2-chloro-4-nitrophenyl-α-D-maltotrioside (CNP-G3) to release 2-chloro-4-nitrophenol and form 2-chloro-4-nitrophenyl-α-D-maltoside (CNP-G2), maltotriose (G3) and glucose (G). The rate of formation of the 2-chloro-4-nitrophenol can be detected spectrophotometrically at 405 nm to give a measurement of α-amylase activity in the sample.

KIT COMPONENTS

For in vitro diagnostic use only.

The components of the kit are stable until expiration date on the label.

Keep away from direct light sources.

**AMY R1 2H120: 6 x 20 ml (liquid) white cap
6U168: 3 x 56 ml (liquid) white cap**

Composition: CNP-G3 2.3 mM, NaCl 350 mM, calcium acetate 6 mM, potassium thiocyanate 600 mM, Good's buffer pH 6.0 100 mM, stabilizers and non-reactive components.

Store all components at 2-8°C.

REAGENT PREPARATION

Use reagent ready to use.

Stability: up to expiration date on labels at 2-8°C.

Stability since first opening of vials: preferably within 60 days at 2-8°C.

PRECAUTIONS

AMY R1: It is not classified as hazardous.

Follow required safety procedures when handling all laboratory reagents.

SPECIMEN

Serum, plasma (heparinate only). Urine.

Amylase is stable in serum and plasma sample up to 2 months at 2-8°C.

EXPECTED VALUES

Serum - plasma: < 96 U/l (< 1.60 µkat/l)
Urine: < 480 U/l (< 8.00 µkat/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA - MULTINORM CHEMA

with normal or close to normal control values

QUANTIPATH CHEMA - MULTIPATH CHEMA

with pathological control values.

If required, a multiparametric, human based calibrator is available:

AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity

the method is linear up to 2000 U/l.

If a ΔA/min of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)

the limit of detection is 0.91 U/l.

Interferences

no interference was observed by the presence of:

hemoglobin ≤ 500 mg/dl
bilirubin ≤ 50 mg/dl
lipids ≤ 1200 mg/dl

Precision

intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	67.89	0.97	1.42
sample 2	171.67	2.61	1.52

inter-assay (n=21)	mean (U/l)	SD (U/l)	CV%
sample 1	67.81	1.93	2.85
sample 2	175.16	4.92	2.81

Methods comparison

a comparison between Chema and a commercially available product gave the following results:

Amylase Chema = x
Amylase competitor = y
n = 155

y = 1.071x - 0.54 U/l r² = 0.997

WASTE DISPOSAL

This product is made to be used in professional laboratories.

P501: Dispose of contents according to national/international regulations.