APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	CL
APP. CODE:	337
WAVELENGTH (Sec/Pri):	700 - 480
ASSAY:	1-POINT      TIME: 10 POINT: 31
SAMPLE VOL:	NORMAL: 2 DECREASE: 2 INCREASE: 5
	R1 VOLUME: 350 R2 VOLUME: 0 R3 VOLUME: 0 R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.0
EXPECTED VALUE:	90 - 110
UNIT:	mEq/l
INSTR. FACTOR (y=ax+b):	a=1    b=0

 Chema Diagnostica  
 Via Campania 2/4  
 60030 Monsano (AN) - ITALY - EU  
 phone +39 0731 605064  
 fax +39 0731 605672  
 e-mail: mail@chema.com  
 website: http://www.chema.com

ITALIANO rev. 16/07/2024

CLORO	
CL 2H500	10 x 50 ml

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**INTENDED USE**  
Reagent for quantitative in vitro determination of chloride in biological fluids.

**PRINCIPLE OF THE METHOD**  
Chloride ions react with mercuric ions, giving available an equal quantity of thiocyanate ions. Thiocyanate ions react with trivalent ferric ions present in solution to form a red colored complex with an absorbance peak at 480 nm.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**CL R1 10 x 50 ml (liquid) white cap**

Composition: mercury(II) thiocyanate 2.2 mM, mercury(II) chloride 0.7 mM, iron (III) nitrate 19 mM.

Store all components at 15-25°C.

**REAGENT PREPARATION**  
Use reagent ready to use.  
Stability: up to expiration date on labels at 15-25°C.  
Stability since first opening of vials: preferably within 60 days at 15-25°C.

**PRECAUTIONS**  
**CL R1:** It is not classified as hazardous.

**SPECIMEN**  
Serum, plasma heparinate. Separation of cells from plasma should be prompt. Sweat is a suitable sample. Use 24 hours urine.  
Dilute sample urine 1:2 with redistilled water and multiply results by two.

**EXPECTED VALUES**  
serum/plasma: 98 - 110 meq/l  
urine: 110 - 250 meq/24h  
(dietary variations are possible) up to 30 meq/l  
sweat: up to 30 meq/l

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**  
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:  
**QUANTINORM CHEMA-MULTINORM CHEMA** with normal or close to normal control values  
**QUANTIPATH CHEMA-MULTIPATH CHEMA** with pathological control values.  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

**AUTOCAL H**  
Please contact Customer Care for further information.

Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Linearity**  
the method is linear up to 200 meq/l.  
If the limit value is exceeded, it is suggested to dilute sample 1+9 with distilled water and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 1.5 meq/l.

**Interferences**  
no interference was observed by the presence of:  
hemoglobin ≤ 500 mg/dl  
bilirubin ≤ 32 mg/dl  
lipids ≤ 500 mg/dl

**Precision**  
intra-assay (n=10) mean (meq/l) SD (meq/l) CV%  
sample 1 114.80 1.48 1.30  
sample 2 111.00 1.41 1.30

inter-assay (n=20) mean (meq/l) SD (meq/l) CV%  
sample 1 117.03 2.95 2.50  
sample 2 113.44 3.26 2.90

**Methods comparison**  
a comparison between Chema and a commercially available product gave the following results:

Chloride Chema = x  
Chloride competitor = y  
n = 83

y = 0.869x + 14.402 meq/l    r<sup>2</sup> = 0.927

**WASTE DISPOSAL**

This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

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