APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912			
TEST:	EPS		
APP. CODE:	333		
WAVELENGTH (Sec/Pri):	700 - 415		
ASSAY:	RATE-A	TIME: 10	POINT: 14 - 25
		DILUENT: water	
SAMPLE VOL:	NORMAL: 10 DECREASE: 8 INCREASE: 12		
	R1 VOLUME: 240 R2 VOLUME: 60 R3 VOLUME: 0 R4 VOLUME: 0	DILUENT: 5	DILUENT: 5
ABS LIMIT:	32000 - INC		
PROZONE LIMIT:	0 - UPPER		
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)		
SD LIMIT:	0.250		
DUPLICATE LIMIT:	3%		
ST. 1 CONC:	0.00		
EXPECTED VALUE:	0 - 100		
UNIT:	U/l		
INSTR. FACTOR (y=ax+b):	a=1	b=0	

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 858)			
TEST NAME:	EPS		
SAMPLE:	Volume 10 µl	Dilution 0 µl	
REAGENTS:	R1 Volume 240 µl R2 Volume 60 µl	Dilution 0 µl Dilution 0 µl	
WAVELENGTH:	Pri. 410 Sec. 700		
METHOD:	RATE		
REACTION SLOPE:	+		
MEASURING POINT 1:	First 17	Last 27	
MEASURING POINT 2:	First	Last	
REAGENT OD LIMIT:	First L -0.1 Last L -0.1	First H 0.5 Last H 0.5	
DYNAMIC RANGE:	L 6	H 1500	
CORRELATION FACTOR:	A 1	B 0	
LINEARITY LIMIT:	15%		
UNIT:	U/l		
CALIBRATION TYPE:	AB		
FORMULA:	Y = AX + B		

 Chema Diagnostica
 Via Campania 2/4
 60030 Monsano (AN) - ITALY - EU
 phone +39 0731 605064
 fax +39 0731 605672
 e-mail: mail@chema.com
 website: http://www.chema.com

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INTENDED USE
 Reagent for quantitative in vitro determination of amylase in biological fluids.

PRINCIPLE OF THE METHOD
 The enzyme α-amylase (EC 3.2.1.1, 1,4 α-D-glucose glucanohydrolase) hydrolyzes the EPS to release several different fragments. The fragments so formed are completely hydrolyzed to 4-nitrophenol and glucose by α-glucosidase. The 4-nitrophenol formed is detected spectrophotometrically at 405 nm to give a measurement of α-amylase activity in the sample.
 The present method has been made according to IFCC.

KIT COMPONENTS
For in vitro diagnostic use only.
 The components of the kit are stable until expiration date on the label.
 Keep away from direct light sources.

AMY EPS R1 2H100: 4 x 20 ml (liquid) white cap
6U280: 4 x 56 ml (liquid) white cap

AMY EPS R2 EA 2H100: 2 x 10 ml (liquid) red cap
EA 6U280: 4 x 14 ml (liquid) red cap

Composition in the test: Hepes buffer 50 mM pH 7.10, NaCl 70 mM, calcium acetate 1.0 mM, EPS-G7 5.0 mM, α-glucosidase 6 KU/l.

Store all components at 2-8°C.

REAGENT PREPARATION
 Use reagent ready to use.
 Stability: up to expiration date on labels at 2-8°C.
 Stability since first opening of vials: preferably within 60 days at 2-8°C.

PRECAUTIONS
AMY EPS R1: Warning. May cause an allergic skin reaction (H317).
 Wear protective gloves (P280) Avoid breathing dust / fume / gas / mist / vapours / spray (P261). If skin irritation or rash occurs: get medical advice / attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364).

AMY EPS R2: Warning. May cause an allergic skin reaction (H317).
 Wear protective gloves (P280). Avoid breathing dust / fume / gas / mist / vapours / spray (P261). If skin irritation or rash occurs: get medical advice / attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364).

SPECIMEN
 Serum, plasma (heparinate only). Urine.
 Amylase is stable in serum and plasma sample up to 2 months at 2-8°C.

EXPECTED VALUES		
Serum - plasma:	28 - 100 U/l	(0.47 - 1.67 µkat/l)
Random urine:	≤ 460 U/l	(≤ 7.68 µkat/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION
 It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA - MULTINORM CHEMA
 with normal or close to normal control values
QUANTIPATH CHEMA - MULTIPATH CHEMA
 with pathological control values.
 If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE
Linearity
 the method is linear up to 1500 U/l.
 If a ΔA/min of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
 the limit of detection is 6 U/l.

Interferences
 no interference was observed by the presence of:
 hemoglobin ≤ 200 mg/dl
 bilirubin ≤ 48 mg/dl
 lipids interfere in low values

Precision			
intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	77.90	0.74	0.90
sample 2	194.80	1.99	1.00

inter-assay (n=20)	mean (U/l)	SD (U/l)	CV%
sample 1	75.77	1.90	2.50
sample 2	194.15	2.39	1.20

Methods comparison
 a comparison between Chema Amylase FL and a commercially available product gave the following results:

$$\begin{aligned}
 & \text{Amylase Chema} = x \\
 & \text{Amylase competitor} = y \\
 & n = 108 \\
 & y = 1.067x + 5.21 \text{ U/l} \quad r^2 = 0.99
 \end{aligned}$$

WASTE DISPOSAL
 This product is made to be used in professional laboratories.
 P501: Dispose of contents according to national/international regulations.

AMYLASE EPS FL	
EA 2H100	4 x 20 + 2 x 10 ml
EA 6U280	4 x 56 + 4 x 14 ml

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